

columns_explained

FIELD	Description	Format	
Recording timestamp	The timestamps that unique identify a recording and which can, later, be used to analyse data in Tobii's proprietary software.	Milliseconds	
Computer timestamp	The Recording timestamp in the computer clock	Milliseconds	
Sensor	Sensor type	Text	
Project name	Control group experiment	Text	
Export date	Date when the Recording was exported	YYYY-MM- DD	
Participant name	Participant name	Text	
Recording name	Recording name	Text	
Recording date	Date when the Recording was performed in this time zone.	YYYY-MM- DD	
Recording date UTC	Date when the Recording was performed in UTC.	YYYY-MM- DD	
Recording start time	Start time of the Recording in this time zone.	YYYY-MM- DD	
Recording start time UTC	Start time of the Recording in UTC format	YYYY-MM- DD	
Recording duration	Total duration of the record	Milliseconds	
Timeline name	Name of the Timeline used during the Recording.	Text	
Recording Fixation filter name	The name of the Fixation Filter applied to the Recording eye tracking data in the export.	Text	
Recording software version	The version of the software used to make the Recording.	Text	
Recording resolution height	Vertical screen resolution used during the Recording.	Pixels	
Recording resolution width	Vertical screen resolution used during the Recording.	Pixels	
Recording monitor latency	The monitor latency setting for the Recording. Stimulus start and end Event timestamps have been off- set by this number to account for the monitor latency.	Milliseconds	
Eyetracker timestamp	The Recording timestamp in the eye tracker clock.	Milliseconds	
Event	Name of the event	Text	
Event value	The event value	Text	
Gaze point X	Horizontal raw gaze coordinates for both eyes combined.	Pixels (MCS)	
Gaze point Y	Vertical raw gaze coordinates for both eyes combined.	Pixels (MCS)	
Gaze point left X	Left eye horizontal raw gaze coordinates.	Pixels (DACS)	
Gaze point left Y	Left eye vertical raw gaze coordinates.	Pixels (DACS)	
Gaze point right X	Right eye horizontal raw gaze coordinates	Pixels (DACS)	
Gaze point right Y	Right eye vertical raw gaze coordinates.	Pixels (DACS)	
Gaze direction left X	Left eye horizontal unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Gaze direction left Y	Left eye vertical unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Gaze direction left Z	Left eye depth unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Gaze direction right X	Right eye horizontal unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Gaze direction right Y	Right eye vertical unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Gaze direction right Z	Right eye depth unit vector for the direction of the gaze,	Normalized coordinates (DACS)	
Pupil diameter left	Estimated size of left pupil	Millimeters	
Pupil diameter right	Estimated size of left pupil	Millimeters	
Validity left	Indicates if the eyes have been correctly identified	Valid / Invalid	
Validity right	Indicates if the eyes have been correctly identified	Valid / Invalid	
Eye position left X (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Eye position left Y (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Eye position left Z (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Eye position right X (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Eye position right Y (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Eye position right Z (DACSmm)	3D position of the eyes	Millimeters (DACS)	
Gaze point left X (DACSmm)	Left eye horizontal raw gaze coordinates	Millimeters (DACS)	
Gaze point left Y (DACSmm)	Left eye vertical raw gaze coordinates	Millimeters (DACS)	
Gaze point right X (DACSmm)	Right eye horizontal raw gaze coordinates	Millimeters (DACS)	
Gaze point right Y (DACSmm)	Right eye vertical raw gaze coordinates	Millimeters (DACS)	
Gaze point X (MCSnorm)	Horizontal raw gaze coordinates on the Media for both eyes combined.	Normalized coordinates (MCS)	
Gaze point Y (MCSnorm)	Vertical raw gaze coordinates on the Media for both eyes combined.	Normalized coordinates (MCS)	
Gaze point left X (MCSnorm)	Left eye horizontal raw gaze coordinates on the Media.	Normalized coordinates (MCS)	
Gaze point left Y (MCSnorm)	Left eye vertical raw gaze coordinates on the Media.	Normalized coordinates (MCS)	
Gaze point right X (MCSnorm)	Right eye horizontal raw gaze coordinates on the Media.	Normalized coordinates (MCS)	
Gaze point right Y (MCSnorm)	Right eye vertical raw gaze coordinates on the Media.	Normalized coordinates (MCS)	
Presented Stimulus name	The name of the Stimulus being presented to the Par- ticipant.	Text	
Presented Media name	The name of the Media presented to the Participant.	Text	
Presented Media width	Width of the Media as presented on the screen to the Participant, including any scaling set in the Stim- ulus properties.	Pixels	
Presented Media height	Height of the Media as presented on the screen to the Participant, including any scaling set in the Sti	Pixels	
Presented Media position X (DACSpx)	The horizontal position of the Media on the screen. The value represents the positions of the top left corner of the Media in relation to the top left corner of the screen.	Pixels (DACS)	
Presented Media position Y (DACSpx)	The vertical position of the Media on the screen. The value represents the positions of the top left corner of the Media in relation to the top left corner of the screen.	Pixels (DACS)	
Original Media width	The original size of the Media presented to the Participant.	Pixels	
Original Media height	The original size of the Media presented to the Participant.	Pixels	
Eye movement type	Type of eye movement classified by the fixation filter settings applied during the gaze data export.	Fixation / Saccade / Unclassified / EyesNotFound	
Gaze event duration	The duration of the cur- rently active eye move- ment.	Milliseconds	
Eye movement type index	Reprensets the order in which an eye movement was recorded. The index is an auto-increments number starting with 1 for each movement type	Number	
Fixation point X	Horizontal coordinate of the fixation point. This is affected by the setting of the Fixation Filter.	Pixels (DACS)	
Fixation point Y	Vertical coordinate of the fixation point. This is affected by the setting of the Fixation Filter.	Pixels (DACS)	
Fixation point X (MCSnorm)	Horizontal coordinate of the fixation point on the media. This is affected by the setting of the Fixation Filter.	Normalize coordinates (MCS)	
Fixation point Y (MCSnorm)	Vertical coordinate of the fixation point on the media. This is affected by the setting of the Fixation Filter.	Normalize coordinates (MCS)	
Mouse position X	The horizontal position of the mouse.	Pixels (DACS)	
Mouse position Y	The vertical position of the mouse.	Pixels (DACS)	